

Make It Right

Centaur Animal Health pledges to **Make It Right** with its customers: current and former. No marketing cliché of “New and Improved”, rather robust scientific improvements tested with field samples. We want to earn your trust by explaining what has been done to address previous issues, and assure you that when you choose the **FP ELISA II**, you can be assured of choosing the best test available from serial to serial.



FP ELISA II

Timely delivery!
Orders shipped
the **same day**,
when placed by
3pm central time.

Call :
1-800-236-6180

**The components are the same
- the procedure has changed.**

✓ **Make it Right**

With more accurate results, the **FP ELISA II** has returned to separate incubations for the sample and the conjugate, with a distilled water wash between. This eliminates the ‘hook’ or prozone effect responsible for several false negatives in 2013, and greatly decreases the number of false positives.

✓ **Make it Right**

By increasing the concentration of conjugate by 20% to make positives more visible.

✓ **Make it Right**

By evaluating each production serial with an expanded panel of problematic samples collected from the field. With every serial, a new filtration step assures no contaminants are introduced.

✓ **Make it Right**

Providing a stop solution that allows you the option of halting color development and allowing interpretation of results at a later time.

✓ **Make it Right**

Former customers will be provided with a **FREE** first kit to evaluate the protocol change and prove to you that we have done what we said we would do.

✓ **Make it Right**

Same pricing we have used for the past 5 years! And with new Competitive Loyalty pricing for labs purchasing more than 20 kits per year, with **NO CONTRACT REQUIRED**.

✓ **Make it Right**

With timely delivery of orders being shipped the same day, when placed by 3 pm central time.

✓ **Make it Right**

The only EIA test able to detect antibodies to the envelope protein of EIAV, as well as the viral core protein.

✓ **Make it Right**

The easiest to interpret test - color in the well equals positive, while absence of color equals negative

The FP ELISA II is manufactured exclusively for Centaur Animal Health by Safepath Laboratories, a USDA licensed establishment.

centaur
ANIMAL HEALTH

1351 West Highway 56 • Olathe, Kansas 66061
913-390-6184 Office • 800-236-6180 Toll Free
centauranimalhealth.com