

# Foalchek<sup>®</sup>

A foalside test for Failure of Passive Transfer

The little one deserves a chance to grow healthy.

Assure it has adequate maternal immunity with Foalchek<sup>®</sup>.

- Available in 10 and 24 test packaging.
- 12-18 months shelf-life allows use across 2 foaling seasons.
- Utilizes serum, plasma, or whole blood.
- Results in less than 6 minutes allows treatment on the same farm visit.
- Test range 200-800 mg/dl
- 81-97% correlation with RID



Available from select distributors.

For more information,  
please contact:

**Centaur Animal Health**

Toll-Free: 800-236-6180

Fax: 913-390-5907

Email:

[sales@centauranimalhealth.com](mailto:sales@centauranimalhealth.com)

[www.centauranimalhealth.com](http://www.centauranimalhealth.com)

**centaur**  
ANIMAL HEALTH

1351-F W. Highway 56 • Olathe, KS 66061