

Canine Heartworm Antigen Test Kit®

**American Heartworm Society “recommends annual
antigen and microfilaria testing”**

(Please see the guidelines at www.heartwormsociety.org)

- Test just **1** drop of serum, plasma or **whole blood**.
- Buffer helps the sample migrate across reaction surface.
- Do not have to stand over the test to activate – can set timer and come back to read results in 10 minutes.
- Available in 5, 15 and 50 test packaging with 12 month shelf-life.
- Accuracy: Sensitivity – 86% Specificity – 96%

For more information contact:

Centaur Animal Health

Toll-Free: 800-236-6180

Fax: 913-390-5907

Email:

sales@centauranimalhealth.com

www.centauranimalhealth.com



Made in USA
By SafePath
for Centaur
Animal Health



centaur
ANIMAL HEALTH

1351-F W. Highway 56 • Olathe, KS 66061