

Canine Parvovirus Antigen Test Kit

- **Collect sample with enclosed swab**
- **Dilution Buffer (included) helps the sample migrate across reaction surface.**
- **Do not have to stand over the test to activate - can set timer and come back to read results in 10 minutes.**
- **Available in 5 and 30 test packaging with 12 month shelf-life.**
- **Accuracy: Sensitivity - 100% Specificity - 99%.**



Made in USA
By SafePath[®]
for Centaur
Animal Health

For more information contact:

centaur
ANIMAL HEALTH

1351-F W. Highway 56 • Olathe, KS 66061
Toll-Free: 800-236-6180 • Fax: 913-390-5907
Email: sales@centauranimalhealth.com • www.centauranimalhealth.com